



Kaye ValProbe Transport/Storage Tray

The ValProbe Transport/Storage Tray is a new accessory **designed to simplify the carrying, storage and management** of the Kaye ValProbe loggers during use.

The ValProbe Transport/Storage Tray can **accommodate up to 20 Kaye ValProbe Loggers** of any type.

Key Features

- High quality, robust design
- Can store up to 20 loggers
- Can be used to securely store the loggers on the shelf
- Good for transportation of loggers after programming
- If you need to program many loggers – it is helpful for sorting non-programmed and pre-programmed loggers to avoid mistakes
- Made with Stainless steel and other materials acceptable for clean room environment, including autoclaves
- Good for Pharma environment, easy to clean with alcohol and other disinfectant sprays
- Rubber feet provides stability when placed on a table or cart
- Size: 300 x 210 x 95 mm (approximate)
- Includes a unique serial number for easy identification



Part Number

V2569

Please contact your local sales representative.

Ask for a quote at: kaye@amphenol-sensors.com

Contact us at: www.kayeinstruments.com/en/contact

www.kayeinstruments.com

© 2017 Amphenol Corporation. All Rights Reserved. Specifications are subject to change without notice. Other company names and product names used in this document are the registered trademarks or trademarks of their respective owners.